

Capstone Design: Preparation of the Subtalar Joint in TTC Arthrodesis

Mackenzie Campbell, Michael De Biasio, Kevin Wang,
Irina Zhu

December 2, 2020



Background: Surgical Procedure

Tibiotalocalcaneal
(TTC) arthrodesis
= fusion




~100mm
incision

=




Joint preparation

+




Ankle Fixation

Background: Subtalar Joint Anatomy



A way to prepare the **subtalar joint** through a **plantar incision** in patients undergoing TTC nailing in a **trauma setting** that improves **bone fusion rates**



A way to prepare the **subtalar joint** through a **plantar incision** in patients undergoing TTC nailing in a **trauma setting** that improves **bone fusion rates**

Space &
Visualization
Constraints

Cost to
Produce
<\$100 CAD




Maximize
surgeon
confidence

Cartilage
and bone
disruption



Minimize
surgical time
& # of steps

Compatible
with bone
graft use


Proposed Solution Overview




Proposed Solution: Guide Tube



Cutting Tool




6mm




10mm max

3mm



100 mm

Cutting Tool



Proposed Solution: Materials and Cost

Materials

Martensitic 400-series stainless steel with blade made of **tungsten carbide**.

Costs


1mm-thick 400-series stainless steel: **USD \$0.815/lb.** (current mkt value).

Tungsten carbide: **USD \$9-12.00/lb.** (approx. mkt value).

Manufacturing: ~ **USD \$30.00/tool.**

Total cost: ~ **CAD \$40.60/tool**

Low fidelity Prototyping






Force Calculations

Friction on hinges	$6.18 \times 10^{-9} \text{ N}$
--------------------	---

Max blade force required to cut cartilage	27.6 N
---	---------------

Handle force	9.17 N
Force on pin	3.93 N

Force on pin (N)



Experimental setup for measuring cartilage cutting force


Blade Selection



Open
curette

Closed
curette

Anterior
lesion
curette



Modified anterior
lesion curette

Force (N)	15		50	
Length to fixed end (mm)	5	10	5	10
Max Stress (N/m ²)	6×10^7		2×10^8	
Max Displacement (μm)	1.7	5.7	2.7	9.1


Tungsten carbide yield strength: $3.35 - 5.30 \times 10^8 \text{ N/m}^2$

Factor of Safety: 2.5-8

Surface Area Calculations

Underside of talus

Top of calcaneus



1

Talus Area: 600 mm²

Calcaneus Area: **726 mm²**


2

With scrapes of size
2mm x 10mm x 1mm:

SA scraped (%)	# of scrapes required
25%	7.5
50%	16.5

~5 minutes

Summary, Design Successes and Failures



Next Steps

- Make a higher-fidelity prototype
 - Give to surgeons, test forces and cartilage scraping
- Gather more feedback from surgeons
- Evaluate subtalar fusion rates

Acknowledgements

Dr. Spencer Montgomery, MD

Dr. Jeremy LaMothe, MD

Dr. Amit Atrey, MD



Dr. Mansur Halai, MD

Dr. Amir Khoshbin, MD

Professor Chris Bouwmeester

Karly Franz

Gary Hoang



References

- K. Yoshimoto, et al., “Does Preparation of the Subtalar Joint for Primary Union Affect Clinical Outcome in Patients Undergoing Intramedullary Nail for Rheumatoid Arthritis of the Hindfoot and Ankle?,” The Journal of Foot & Ankle Surgery, vol. 59, no. 5, pp. 984-987, 2020. Available: [https://www.jfas.org/article/S1067-2516\(20\)30165-4/fulltext](https://www.jfas.org/article/S1067-2516(20)30165-4/fulltext)
- C. L. Brockett and G. J. Chapman, “Biomechanics of the ankle,” Orthopaedics and Trauma, vol. 30, no. 3, pp. 232-238, 2016. Available: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4994968/>.
- G. E. Quill, “Tibiototalcalcaneal Arthrodesis With Medullary Rod Fixation,” Techniques in Foot and Ankle Surgery, vol. 2, no. 2, pp. 135-143, 2003. Available: https://journals.lww.com/techfootankle/Abstract/2003/06000/Tibiototalcalcaneal_Arthrodesis_With_Medullary_Rod.9.aspx.
- B. H. Lee, C. Fang, R. Kunnasegaran, and G. Thevendran, “Tibiototalcalcaneal Arthrodesis With the Hindfoot Arthrodesis Nail: A Prospective Consecutive Series From a Single Institution,” The Journal of Foot & Ankle Surgery, vol. 57, no. 1, pp. 23-30, 2018. Available: [https://www.jfas.org/article/S1067-2516\(17\)30358-7/fulltext](https://www.jfas.org/article/S1067-2516(17)30358-7/fulltext).
- K. Akiyama, T. Sakai, N. Sugimoto, H. Yoshikawa, and K. Sugamoto, “Three-dimensional distribution of articular cartilage thickness in the elderly talus and calcaneus analyzing the subchondral bone plate density,” Osteoarthritis and Cartilage, vol. 20, no. 4, pp. 296-304, 2012. Available: <http://www.sciencedirect.com/science/article/pii/S1063458412000209>.
- Stryker, “Ankle Arthrodesis Nail: Operative Technique”, pp. 1-31, 2009. Available: <https://www.strykermeded.com/media/1602/t2-ankle-arthrodesis-nail.pdf>.



References

- P. Chen, J. Sui, and C. Wang, “Cutting Force Analysis of Bovine Acetabular Cartilage,” *Procedia CIRP*, vol. 89, pp. 189-193, 2020. Available: <https://www.sciencedirect.com/science/article/pii/S2212827120305163>.
- M. Mahvash, L.M. Voo, D. Kim, K. Jeung, J. Wainer, and A.M. Okamura, “Modeling the forces of cutting with scissors,” *IEEE transactions on bio-medical engineering*, vol. 55, no. 3, pp. 848-56, 2008. Available: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2709828/>.
- G. L. Bennett, B. Cameron, G. Njus, M. Saunders, and D. Kay, “Tibiototalcaneal Arthrodesis: A Biomechanical Assessment of Stability,” *Foot & Ankle International*, vol. 26, no. 7, pp. 530-536, 2005. Available: <https://pubmed.ncbi.nlm.nih.gov/16045843/>.
- Y. Yasui, C. P. Hannon, D. Seow, and J. G. Kennedy, “Ankle arthrodesis: A systematic approach and review of the literature,” *World Journal of Orthopaedics*, vol. 7, no. 11, pp. 700–708, 2016. Available: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5112338/>.
- M.M. Van Meter and R.A. Adam, “Costs associated with instrument sterilization in gynecologic surgery”, *American Journal of Obstetrics & Gynecology*, vol. 215, no. 5, 2016. Available: [https://www.ajog.org/article/S0002-9378\(16\)30358-1/pdf](https://www.ajog.org/article/S0002-9378(16)30358-1/pdf).
- “ACUFEX Currettes,” Smith & Nephew. Available: <https://www.smith-nephew.com/canada/products/extremities-and-limb-restoration/acufex-currettes/>

